

# SPE-CapNMR™ Application Note


By: Aaron Wilson and David Detlefsen

A basic SPE process is described that can enrich and concentrate a sample in deuterated solvent for analysis by the CapNMR™ Probe.


In the first case, Amitriptyline is used as a model compound in a basic SPE scheme to enrich the concentration of the sample. In the second case, Buspirone is used in a recovery scheme.

# Sample Trapping Conditions



- Amitriptyline (MW = 314 g/mole)
  - Use Michrom BioResources microtrap
  - Pre-wash with 250  $\mu\text{L}$   $\text{H}_2\text{O}$
  - Condition with 500  $\mu\text{L}$  Methanol
  - Equilibrate with 500  $\mu\text{L}$  of  $\text{D}_2\text{O}$  with 1.2% HFBA
  - Load 500  $\mu\text{L}$  of Amitriptyline (20  $\text{ng}/\mu\text{L}$ ; 64  $\mu\text{M}$ ) in  $\text{H}_2\text{O}$
  - Elute three 8- $\mu\text{L}$  fractions of Methanol- $\text{d}_4$
  - CapNMR<sup>TM</sup> Analysis
    - Triplet at 5.7 PPM



- Buspirone (MW = 386 g/mole)
  - Use Michrom BioResources microTrap
  - Clean trap with 100  $\mu\text{L}$  acetonitrile
  - Equilibrate with 100  $\mu\text{L}$  sample solvent with 0.05% HFBA
  - Load sample dissolve in solvent with 0.05% HFBA
  - Dry with 500  $\mu\text{L}$  air injections
  - Elute with 15  $\mu\text{L}$  Acetonitrile-  $\text{d}_3$
  - CapNMR<sup>TM</sup> Analysis
    - Doublet at 8.39



# Sample Trapping - Amitriptyline






## Before Trapping

- Amitriptyline (256 Scans in 8.5 min)
- 20 ng/ $\mu$ L sample solution
- 8  $\mu$ L manually injected into probe
- $\sim 0.1$   $\mu$ g in  $V_{\text{obs}}$  of CapNMR<sup>TM</sup> probe
- S/N = 6.5

## After Trapping

- Amitriptyline (256 Scans in 8.5 min)
- Other conditions same as above
- SPE fraction: 2<sup>nd</sup> 8- $\mu$ L elution
- 8  $\mu$ L manually injected into probe
- $\sim 2$   $\mu$ g in  $V_{\text{obs}}$  of CapNMR<sup>TM</sup> probe
- S/N = 132 (20 fold enrichment)

# Sample Trapping - Buspirone



- Reference: Trap from 90:10
  - 20 ug loaded onto trap
  - 10  $\mu$ L of 5.2 mM solution
  - 128 scan in 4.5 min at 600 MHz
  - 10  $\mu$ L manually injected into probe
  - Relative peak area 9.86 (100%)
- Trial 1: Trap from 90:10
  - 5 ug loaded onto trap
  - 250 ul of 52 uM solution
  - 128 scan in 4.5 min at 600 MHz
  - 10  $\mu$ L manually injected into probe
  - Relative peak area 1.04 (42% of reference)
- Trial 2: Trap from 98:2
  - 5 ug loaded onto trap
  - 250  $\mu$ L of 52 uM solution
  - 128 scan in 4.5 min at 600 MHz
  - 10  $\mu$ L manually injected into probe
  - Relative peak area 2.95 (120% of reference)

# Sample Trapping Summary

- Preconcentration simplifies trace NMR analysis
  - Presents opportunities for on-line preconcentration
  - Michrom Trap Cartridges
    - enable microscale sample preconcentration
- Efficient Sample Preconcentration
  - 20 fold enrichment
  - 50% or greater recovery
- Trapping Solution Dependencies
  - Aqueous:organic ratio affects recovery
  - HFBA ion pairing agent important