

**TXI Probe Characterization Sample**  
**14 Jan 2008; DX**

NMR Characterization Reagents for TXI CapNMR Probes

Vendor: Cambridge Isotope Labs

(1) Methanol- $p_4$  ( $^{13}\text{C}$ , 99%); Catalog #CLM-359

- 1 g is \$125
- For a 200 mM stock solution in DMSO- $d_6$ , a 5-mL volume contains 1 mmole of  $^{13}\text{C}$ -methanol- $p_4$  (33 mg or 40.5  $\mu\text{L}$  @ 25° C. This roughly costs \$4 (plus the cost of the solvent).
- Seal remaining methanol- $p_4$  in a small vial using special black cap, and Parafilm; keep dry.

(2) Urea ( $^{15}\text{N}_2$ , 98%); Catalog #NLM-233

- 1 g is \$166
- For a 200 mM stock solution in DMSO- $d_6$ , a 5-mL volume contains 1 mmole of  $^{15}\text{N}$ -urea (62 mg). This roughly costs \$10.30 (plus the cost of the solvent).
- Keep powder dry.

The total cost of the 5-mL solution is less than \$20.

### Formula Summary

- Measure 41  $\mu\text{L}$  of labeled methanol (using syringe) and add to 5-mL volumetric flask
- Weigh 62 mg (0.0620 g) of labeled urea (using analytical balance) and add to 5-mL volumetric flask
- Dissolve both to 5.00 mL of DMSO- $d_6$  in the 5-mL volumetric flask

Make 10 aliquots of about 0.5 mL each into small glass vials sealed with the “special” black caps and Parafilm.

The final solution is dissolved in DMSO- $d_6$  and contains (asterisks indicate isotopic label location):

- 200 mM  $^{13}\text{C}$ -labeled methanol- $p_4$  ( $\text{C}^*\text{H}_3\text{OH}$ )
- 200 mM  $^{15}\text{N}$ -labeled urea, which is  $\text{H}_2\text{N}^*-(\text{C}=\text{O})-\text{N}^*\text{H}_2$

### Notes

- DMSO- $d_6$  is chosen as a solvent to disallow loss of the amino protons by exchange with the solvent.
- Store labeled reagents well-sealed in the refrigerator.