

Technical Bulletin #C0060: gHMBC with CapNMR

This bulletin describes how to use the CapNMR probe for Ghmhc (^1H - ^{13}C) on a Varian console/system.

Begin by loading an appropriate 'Proton' sequence. Make sure to use the appropriate power levels (see Figure 1). The solenoid microcoil in the CapNMR probe requires less transmitter power than traditional probes. The standard pulse width (pw90) is calibrated at installation for a transmitter power of 45 dB (tpwr=45).

Parameter	Varian Range	Bruker Range
Lock power	1 to 8	-55 to -45
Lock gain	34 to 48	110 to 130
Transmitter power level (dB)	40 to 45	15 to 20
90° Transmitter pulse width (µsec)	2.5 to 5	2.5 to 5
90° Decoupling pulse width (µsec)	5 to 12	7 to 12

Figure 1. Appropriate power levels for the CapNMR Probe

Once there is a good Proton spectrum, type the macro 'Ghmhc' into the command prompt and this will load a slightly modified gHMBC pulse sequence. You can view the pulse sequence with a dps command. The sequence should look like the sequence in Figure 2.




Figure 2. Ghmhc Pulse Sequence

However, the pulse width is probe/console specific so your pulse width values (pw & pwx) should be different than the 2.6 and 4.8 values shown in Figure 2 and Figure 3. The pw90 for proton was calibrated with an array of the pw parameter in a Proton sequence. For tpwr=45, the pw90 for this probe was 2.6 microseconds. The pwx90 for carbon was calibrated with a pwxcal pulse sequence and for pwxlvl=45, the pwx90 for this probe was 4.8 microseconds.

Here is an example of the parameter set used for the Ghmbc pulse sequence. In order to view this window, click on Process and select Text Output and then type 'dg'.

	ACQUISITION	TRANSMITTER	HMBC	PROCESSING
Default	sw 5584.6	tn H1	jlxh 140.0	sb 0.064
Weighting	at 0.128	sfrq 599.763	jnxh 8.0	sbs not used
Display	np 1430	tof -311.8	GRADIENTS	fn 4096
More 2D	bs 32	tpwr 45	gzlvl1 888	2D PROCESSING
Integration	ss 32	pw 2.650	gt1 0.001000	sb1 0.006
Cursors/Line Lists	dl 1.000	DECOUPLER	gzlvl3 444	sbs1 not used
Plot	nt 88	dn C13	gt3 0.001000	procl lp
Text Output	ct 88	dof 2339.0	gstab 0.000500	fn1 2048
	2D ACQUISITION	pwxlvl 45	hsglvl 888	SAMPLE
	swl 36199.1	pwx 4.800	hsgt 0.002000	date Sep 14 2006
	ni 200	dm nnn	SPECIAL	solvent CDC13
	PRESATURATION	dmf 22000	temp not used	sample
	satmode nnn	dpwr 23	spin not used	
	satdly 0	gain	54	

Figure 3. Ghmbc Parameters for the CapNMR Probe

IMPORTANT: GRADIENT POWER LEVELS

The CapNMR probe has a much stronger z-gradient than most standard tube probes. Most common tube probes are 50-80 gauss/cm at 10 amps while the CapNMR is commonly 350-400 gauss/cm at 10 amps. The default gzlvl1 and hsglvl values from Varian are 2000. You may have seen these values as high as 4500 for some tube probes. However, due to a more efficient gradient coil, not as much gradient power (gzlvl) is required for the same gradient field.

A standard gradient value (gzlvl1 & hsglvl) for the CapNMR probe is 700 to 900.

A standard gradient value (gzlvl3) for the CapNMR probe is 350 to 450.

These gradient and homospoil gradient values are automatically calculated from the gcal value in the probe file. The gcal value is calibrated at installation with a profile sequence, which performs a z gradient echo, and is saved in the probe file. Make sure your probe is set correctly (check with 'probe?') and that the 'probegcal' value is correct.

Sample and Data Set:

The sample used was 10% (v/v) 2-ethyl-1-indanone in CDCl₃. Pure 2-ethyl-1-indanone is a clear, brown liquid with a T1 relaxation of approximately 0.8 seconds. This data set was taken on a **5- μ L of 10% 2-ethyl-1-indanone (519 μ g; 3.2 μ mol; v/v in CDCl₃)**, which was injected and then pushed into a 5 microliter flowcell with more CDCl₃ solvent. This data was taken on a Varian Inova at 600 MHz with a TXI (Triple Inverse ¹H, ²H, ¹³C, ¹⁵N) CapNMR probe with a 5 microliter enhanced flowcell in 4 hours.




Figure 4. Ghmbc Spectrum from 10% 2-ethyl-1-indanone




Figure 5. Ghmbc Spectrum (Aliphatic Region) from 10% 2-ethyl-1-indanone




Figure 65. Ghmbc Spectrum (Aromatic Region) from 10% 2-ethyl-1-indanone

Injection Method

A Hamilton Gastight 25- μ L syringe was filled with clean CDCl_3 solvent and used to rinse the probe. Then, 5- μ L of 10% 2-ethyl-1-indanone (519 μ g; 3.2 μ mol; v/v in CDCl_3) was drawn into the syringe and injected into the probe. Using the same syringe, 12 μ L of clean CDCl_3 was picked up and injected into the probe to push the 5- μ L sample into the NMR flow cell of the CapNMR probe. The 12- μ L Push Volume was calibrated in advance of sample injection.

Screenshots

All of the screenshots shown in this Technical Bulletin are from a Varian Inova console running Varian VnmrJ 2.1B software on a Unix workstation. The magnet was an Oxford 600 MHz magnet. (For more information on the Varian VnmrJ software, see <http://www.varianinc.com/cgi-bin/nav?products/nmr/>)